# Be Alive ! DIO UV Active





# **CONTENTS**



# Simple & Powerful Body Design

• Tapered-Straight Combined Shape leads easily placement and increasing initial fixation power and stability

# **UFII Implant Characteristics 1**



n

# **UFII Implant Characteristics 2**

# Hybrid Type SLA Implant

• Hybrid Sand Blast & Acid Etched Surface Treatment Minimize peri-implantitis and stimulate osseointegration





## Minimize Peri-Implantitis

Less friction around cortical bone. Minimizes bone heating

## Stimulate Osseointegration

With ideal Ra value promotes not only rapid osseointegration but also long-term stability



# **Easy & Simple Drilling Protocol** UFII Surgical Kit [UF07]



# **UFII Surgical Drill**

- Easy depth control
- Minimizes bone heating from significant cutting





# **Biological aging**

- 1 The freshly manufactured implant initially has hydrophilic property on its surface.
- 2 When exposed to air, the SLAtreated surface forms a layer of titanium dioxide(TiO<sub>2</sub>), a photocatalyst medium.
- **3** Biological aging occurs as soon as implants are exposed to air.
- 4 Attachment and accumulation of hydrocarbons are on the surface.
- 5 Conversion of surface to hydrophobic and electronegative elements.
- 6 Decreased osseointegration efficiency.



# **UV TREATED UV** Photofunctionalization

- 1 UV photocatalytic activity proceeds when UVC is irradiated on the TiO<sub>2</sub> layer of the implant surface.
- 2 OH-radicals remove hydrocarbons from the implant surface and increases BIC.
- 3 Conversion of a surface from hydrophobic to superhydrophilic induces improved blood wettability.
- 4 Surface modification to electropositive leads to the direct electrostatic interaction of UV treated surface with cells.
- **5** Induces fast bone formation and attains early osseointegration.

Direct Electrostatic Interaction

# **Experimental Results** ①

# **Change Hydrophilicity and Increase Wettability**

- Superhydrophilicity from OH Radical maximizes blood wettability and improves the adhesion of proteins and osteoblasts, speeding up initial bone formation.
- Faster loading time and shorter unstable period after placement.



# Hydrocarbons Matter Removed From Implant Surface

• UV photofunctionalization could remove about half of unavoidably contaminated hydrocarbons on implant surface before placement, which increases BIC ratio.



## **XPS Analysis Result**

By Center for Research Facilities at Pukyung National University, Korea Specimen : Titanium disk (Ø10 \* 3t)



## **Carbon Weight Ratio**



## **Oxigen Weight Ratio**

# **Experimental Results** ③

# **Remove Hydrocarbons and Increase BIC**

- UV treated implant that removes organic matter from the surface induces stable bone formation because the implant has a high BIC (Bone to Implant Contact).
- Implant failure rate is low due to the high level of bone intensity during formation.

# **Experimental Results** ④

# **Earlier Bone Formation**

• UV photofunctionalization converts the implant surface from hydrophobic to superhydrophilic and removes organic matters from the surface, therefore, this modification could accelerate osseointegration and induce early loading of the prosthesis.



# **BIC experiment after UV treatment**

By College of Dentistry at Kyungpook University, Korea





# In Vitro Test Cell Proliferation Experiments After UV Treatment

## By College of Dentistry at Kyungpook University, Korea

- Test specimen: Murine pre-osteoblast MC3T3-E1 cells
- UV treatment induces super hydrophilicity and increases osteoblast proliferation, attachment, and mineralization.

# **UV Activator**

## A 20-second Hydrophilic Surface Treatment Solution

# UV untreatedUV treated1 DayImage: Comparison of the second sec

7 Days

MC3T3-E1 Proliferation

**X 200** Rate of cell growth of UV treated Ti disk increases immensely from day 3.





Bone formation on UV treated Ti disk dramatically increased between 2~3 weeks, compared to SLA's 3~4 weeks

![](_page_6_Picture_13.jpeg)

# **UV Activator**

# Attractive to water, Attractive to patients

- The patented cutting-edge technology
- The world's first UV treatment machine in less than a minute

![](_page_7_Picture_5.jpeg)

## 1 Fast and Powerful

- Induces photocatalytic effect on implant surface with just 20 seconds of UV irradiation.
- The all-new 360° cylindrical form factor and high power increased the absorption rate of UV rays by thoroughly irradiating the titanium surface, achieving superior hydrophilicity dramatically boosting osseointegration.

## 2 Safe and Confident

- Unique double-layer packaging allows UV irradiation without opening an ampoule, for easier usage and to prevent even the slightest chance of contamination
- More then 700,000 UV Active delivered to over 20 countries around the world
- Won iF Design Award 2020, which is one of the world's top three design awards.

![](_page_7_Picture_13.jpeg)

# When do you need UV Active?

# To secure initial stability for weak bone conditions

- When you want to get strong initial stability without bone graft by Sejong Star Dental Clinic, Dr. JongHwan Park
- When bone graft is needed on the wide-area of the maxillary sinus by Yedam Dental Hospital, Dr. JaeSeok Kang
- When the implant is placed immediately after extraction by Michigan Dental Clinic, Dr. Hyangryeon Lee
- When you want to replace the failed implant in the same hole by Yedam Dental Hospital, Dr. JaeSeok Kang

### To complete treatment rapidly 2

- When you want to greatly shorten the treatment period for final restoration by Essence Dental Clinic, Dr. SungBum Hong
- When you have patients who can not visit the clinic frequently by JaeHyun Ahn Dental Clinic, Dr. jaeHyun Ahn

### To differentiate from other clinics 3

- When you want to distinguish your clinic from others by M Plant Dental Clinic, Dr. MinWoo Lee
- When you want to differentiate your implant from other implants (Biological-Aging) is required by Yedam Dental Hospital, Dr. JaeSeok Kang
- When you need to satisfy your patients' needs who want more special implants by M Plant Dental Clinic, Dr. MinWoo Lee

![](_page_7_Picture_27.jpeg)

# UV Active is an effective implant for:

Patients who have soft, weak or insufficient dental bone or need bone graft.

Patients who require immediate implant placement after extraction.

**Patients who require** complex surgeries such as sinus lift or edentulous case patients.

Patients who are too busy to visit the hospital often and desires fast recovery.

> **Elderly patients who** have bad dental bone conditions.

**Patients who require** re-surgery due to implant failure.

# **Testimonials**

![](_page_8_Picture_10.jpeg)

## **Dr. Jaeseok Kang**

### Yedam Dental Hospital, Republic of Korea

Placed UV Active in the wide-area of maxillary sinus and measured ISQ. I saw rapid osseointegration after 2 months.

Strongly recommend UV Active when you need to place a lot of bone graft in the maxillary sinus or if it is just a GBR case.

![](_page_8_Picture_15.jpeg)

## **Dr. Riley Clark**

## WhiteCap Institute, USA

One of my personal game-changers is UV Active which is super-hydrophilic Implant.

We have that special obligation in the doctor-patient relationship to always offer the most contemporary and best practices and although this UV light activation photo functionalization is a newer topic it's an exciting topic and we're excited to dive into it and talk about some of the benefits for us as a clinician.

![](_page_8_Picture_22.jpeg)

## **Dr. Jungwook Seo**

## Yonsei Dental Clinic, Republic of Korea

The benefits of the UV Active are it increases the initial stability when placed after UV surface treatment and it also extends the implant shelf life.

![](_page_8_Picture_26.jpeg)

## Dr. Hyangryeon Lee

## Michigan Dental Clinic, Republic of Korea

UV Active is perfect to use when placing an immediate temporary after extraction or for a GBR case due to its fast osseointegration rate.

My patients are very satisfied not to mention it is also very beneficial for dentists.

![](_page_9_Picture_0.jpeg)

# **Clinical Utilization** of UV Activated Implant

![](_page_9_Picture_2.jpeg)

Case 2 **Clinical Case Report** 

**UV Activated Implant** 

Dr. JaeSeok Kang, Yedam Dental Hospital

Patient Information: Male, 72 years Treatment Plan: #11,13,23 extracted & #11~13, #21~23 Bridge

![](_page_9_Picture_7.jpeg)

![](_page_9_Picture_8.jpeg)

Panoramic view at Initial Examination

Panoramic photo and oral cavity photo after the final prosthesis

Dental formula No.	Implant Size	Bone Density	ISQ					
			ОР	1Week	2Week	3Week	4Week	5Week
#13	DIO HSA Ø3.8 X 10	D2~D3	65	65	75	75	76	75
#11	DIO HSA Ø3.8 X 10	D3	65	66	80	80	81	82
#21	DIO HSA Ø3.8 X 10	D3	68	68	73	73	73	74
#23	DIO HSA Ø3.8 X 10	D2~D3	80	80	82	82	83	83

ISQ measurement value per every week and after surgery

![](_page_9_Picture_13.jpeg)

Post OP

The ISQ trend of implants placed in the anterior maxilla is observed to increase in the 2nd week, with DIOnavi. & UV Active, it was possible to achieve osseointegration in a short period of time.

Treatment Plan: #26~27 Bridge / #27 required extensive alveolar bone-graft at the same time

![](_page_9_Picture_17.jpeg)

Patient Information: Male, 64 years

Pre OP

Post OP

![](_page_9_Picture_20.jpeg)

Dontal formula	ISQ						
No.	ОР	1Week	2Week	3Week			
#26	67	75	81	83			
#27	21	52	64	75			

ISQ measurement value per every week and after surgery

In general, when extensive bone loss has progressed in the maxillary posterior, an implant is placed along with a bone graft and osseointegration time is at least 4 to 5 months. This case took a total of 3 months and the UV-irradiated implant showed a much faster osseointegration rate.

# Overcoming difficult cases by using

![](_page_9_Picture_26.jpeg)

![](_page_9_Picture_27.jpeg)

Final Restoration

4Week
85
81

# Milestones of UV Technologies and Digital Solutions

700K UV Active Implants delivered to more than 20 countries

![](_page_10_Picture_2.jpeg)

# 2010

Modeless Customized Abutment I-FIT (iTero Scanner Base)

# 2011

**UV Active** Surface

# 2012

UFII HSA Implant 3Shape Trios

# 2014

DIOnavi. Digital Surgical Guide System

# 2015 Digital

Abutment

# 2018 **3D** Printing Solutions

UV Activator 1 **UV Active Implant** 

Modeless CAD/CAM Technology Based on Intra Oral Scanner | 12 Years

![](_page_10_Picture_17.jpeg)

2019 UV Activator2 Launched \* Won iF Design Award 2020

![](_page_10_Picture_19.jpeg)

Surface Treatment Launched

2019 **UV Activator 2** DESIGN AWARD 2020

![](_page_10_Picture_23.jpeg)

DIOnavi.Launched

DIOnavi. Full Arch Digital edentulous rehabilitation system

# **Fully Digitalized Implant Solution | 8 Years**

# UV Technology | 11 Years

# Proposal of Implant Treatment Line-up

The strategic line-up offers a simple guide for a patient's decision-making process

![](_page_11_Picture_2.jpeg)

![](_page_12_Picture_1.jpeg)