


UFII Master Kit

Kit Code | UF 06

- Kit composed of drill with outstanding cutting force and durability, etc.
- UFII Narrow $\varnothing 3.0$ $\varnothing 3.3$ / UFII Regular $\varnothing 3.8$ $\varnothing 4.0$ $\varnothing 4.5$ $\varnothing 5.0$ Fixture

Unit mm | Drill Scale 1 : 0.9



UFII Surgical Kit


Kit Code | UF 07

- Kit composed of drill with outstanding cutting force and durability, etc.
- UFII Narrow $\varnothing 3.0$ $\varnothing 3.3$ / UFII Regular $\varnothing 3.8$ $\varnothing 4.0$ $\varnothing 4.5$ $\varnothing 5.0$ $\varnothing 5.5$ Fixture

Unit mm | Drill Scale 1 : 0.9

IMPLANT SYSTEM

UFII INT. SUBMERGED SYSTEM | UFII Master-Surgical Kit



Torque Wrench
DTW 0060



Open Wrench
OW 004



Depth Gauge
DPG 2050


Surgical Tool

Starter Drill

K

Code

SD 1530



Unit mm | Scale 1 : 1.15

Initial Drill


K

D Ø2.0

· Only use UFII Master Kit

Unit mm | Scale 1 : 1

Code	DSU 2007D	DSU 2008D	DSU 2010D	DSU 2011D	DSU 2013D
------	-----------	-----------	-----------	-----------	-----------




K

D Ø2.0

· Only use UFII Surgical Kit

Unit mm | Scale 1 : 1.2

Code	DHF 2015SD
------	------------



Drill stopper

K

D Ø2.0

· Only use UFII Surgical Kit

Unit mm | Scale 1 : 1

Code	ST 2015F	ST 2013F	ST 2011F	ST 2010F	ST 2008F
------	----------	----------	----------	----------	----------








● Sold separately

Parallel Pin

● Sold separately

Unit mm | Scale 1:1


K					
Code	PP 1223	PP 1523	PP 2015	PP 2022	PP 2023
					

Path Pin

● Sold separately

Unit mm | Scale 1:1.2


K	
Code	DPP 3512



Lindemann Drill

Unit mm | Scale 1:1.2

K	D Ø2.0
Code	RMH 2014S



Positioning Guide

● Sold separately

Unit mm | Scale 1:1.2

K	D Ø2.0
Code	PG 0060



Angled Path Pin

● Sold separately

Unit mm | Scale 1:1.15

K		
Code	DAP 4515A	DAP 4515B




Type

Hex A

Hex B

Surgical Tool

Drill



K D Ø2.6

● Sold separately

Unit mm | Scale 1:1

Code	DSU 2707D	DSU 2708D	DSU 2710D	DSU 2711D	DSU 2713D	DSU 2715D	DSU 2718D
Length	7	8.5	10	11.5	13	15	18

K D Ø2.9

Unit mm | Scale 1:1

Code	DSU 3007D	DSU 3008D	DSU 3010D	DSU 3011D	DSU 3013D	DSU 3015D	DSU 3018D
Length	7	8.5	10	11.5	13	15	18

K D Ø3.1

Unit mm | Scale 1:1

Code	DSU3207D	DSU3208D	DSU 3210D	DSU 3211D	DSU 3213D	DSU 3215D	DSU 3218D
Length	7	8.5	10	11.5	13	15	18

K D Ø3.7

Unit mm | Scale 1:1

Code	DSU 3807D	DSU 3808D	DSU 3810D	DSU 3811D	DSU 3813D	DSU 3815D	DSU 3818D
							
Length	7	8.5	10	11.5	13	15	18

K D Ø4.2

Unit mm | Scale 1:1

Code	DSU4307D	DSU4308D	DSU 4310D	DSU 4311D	DSU 4313D	DSU 4315D	DSU 4318D
							
Length	7	8.5	10	11.5	13	15	18

K D Ø4.7

Unit mm | Scale 1:1

Code	DSU 4808D	DSU 4810D	DSU 4811D	DSU 4813D	DSU 4815D	DSU 4818D
						
Length	8.5	10	11.5	13	15	18

Surgical Tool

Profile Drill

K

Unit mm | Scale 1:1

Fixture Size	Ø3.8	Ø4.0	Ø4.5	Ø5.0	Ø5.5
Code	UPD 3805D	UPD 4005D	UPD 4505D	UPD 5005D	UPD 5505D

UPD 3805D: Ø3.8, 16mm length, 30mm total length, blue color.

UPD 4005D: Ø4.0, red color.

UPD 4505D: Ø4.5, grey color.

UPD 5005D: Ø5.0, green color.

UPD 5505D: Ø5.5, white color.

Color

●

●

●

●

○

Tap Drill

K

Unit mm | Scale 1:1

Fixture Size	Ø3.0	Ø3.3	Ø3.8	Ø4.0	Ø4.5	Ø5.0	Ø5.5
Code	UNT D 3015	UNT D 3315	UT (II) 3815	UT (II) 4015	UT (II) 4515	UT (II) 5015	UT (II) 5515

UNT D 3015: Ø3.0, purple color.

UNT D 3315: Ø3.3, yellow color.

UT (II) 3815: Ø3.8, blue color.

UT (II) 4015: Ø4.0, red color.

UT (II) 4515: Ø4.5, grey color.

UT (II) 5015: Ø5.0, green color.

UT (II) 5515: Ø5.5, white color.

Color

●

●

●

●

●

●

○

Drill Extension

K

Unit mm | Scale 1:1.2

DE 3811: Drill Extension, silver color.



Code	DE 3811
------	---------

Fixture Driver Narrow




● Sold separately

K Machine

Unit mm | Scale 1 : 1

Code	UNFDC 3006	UNFDC 3012	UNFDC 3018
			

K Ratchet



Code	UNFDW 3006	UNFDW 3012	UNFDW 3018
			

Fixture Driver Regular / Wide




● Sold separately

K Machine

Unit mm | Scale 1 : 1

Code	FDC 3506	FDC 3512	FDC 3518
			

K Ratchet

Code	FDW 3506A	FDW 3512A	FDW 3518A
			


Surgical Tool

Screw Driver Machine

● Sold separately

K Slot 0.5

Unit mm | Scale 1 : 1

Code	MD 0522	MD 0530
		


K Hex 0.9

Code	MD 0922	MD 0930
		

K Hex 1.2

Code	MD 1219	MD 1222	MD 1230	MD 1234	MD 1239
					

K Torx 1.7



Code	MD 1719	MD 1722	MD 1730
			

Screw Driver Ratchet




● Sold separately

K Slot 0.5







Unit mm | Scale 1 : 1

Code	HD 0510A	HD 0515A
		






K Hex 0.9

Code	HD 0910A	HD 0915A	HD 0920A
			

K Hex 1.2

Code	HD 1205A	HD 1210A	HD 1215A	HD 1220A	HD 1225A	HD 1230A
						

K Torx 1.7

Code	HD 1705A	HD 1710A	HD 1715A	HD 1720A	HD 1725A
					

Surgical Tool


Fixture Mount Driver

· Only use UF II Master Kit

K Hex 2.5		Unit mm Scale 1 : 1	
Type	Machine		Ratchet
Code	MHDC 2520	MHDC 2525	MHDR 2513A MHDR 2518A
			

Solid Abutment Driver

· Only use UF II Master Kit

K		Unit mm Scale 1 : 1.2	
Length	6		12
			
Diameter Ø4.5	HDS 4506A	HDS 4512A	
Ø5.5	HDS 5506A	HDS 5512A	
Ø6.5	HDS 6506A	HDS 6512A	

Ball Abutment Driver

· Only use UF II Master Kit

● Sold separately


K

Unit mm | Scale 1:1

Code

HD 2406A

HD 2412A



Abutment Removing Driver

· Only use UF II Master Kit



K

Unit mm | Scale 1:1

Code


HARD 20315A

HARD 20320A



Surgical Tool

Depth Gauge



K D Ø2.0

Unit mm | Scale 1 : 1

Code

DPG 2050

Open Wrench

- Used to remove Mount if the bone is soft.
- 30 ° Neck angle for convenience in intraoral insertion.



K

Unit mm | Scale 1 : 1

Code

OW 004

Torque Wrench



K

Unit mm | Scale 1 : 1

Code


DTW 0060

UFII Narrow Kit

Kit Code | UF 08


- Kit composed of drill with outstanding cutting force and durability, etc.
- UF II Narrow Ø3.0 Ø3.3 Fixture

Unit mm | Drill Scale 1 : 0.9




Surgical Tool

Starter Drill



K		Unit mm Scale 1 : 1.2	
Code		SD 1530	

Initial Drill





K D Ø2.0		Unit mm Scale 1 : 1.2	
Code		DHF 2015SD	

Drill stopper

K D Ø2.0		Unit mm Scale 1 : 1			
Code	ST 2015F	ST 2013F	ST 2011F	ST 2010F	ST 2008F
					

Parallel Pin

K D Ø2.0		Unit mm Scale 1 : 1	
Code		PP 2022	PP 2023



Surgical Tool


Drill

K

D Ø2.6

Unit mm | Scale 1 : 1

Code	DSU 2708D	DSU 2710D	DSU 2711D	DSU 2713D
------	-----------	-----------	-----------	-----------











Tap Drill

K

Unit mm | Scale 1 : 1

Fixture Size	Ø3.0	Ø3.3
Code	UNTD 3015	UNTD 3315





Fixture Driver

· Maximum Tightening Torque : 50Ncm

 Sold separately

K


Unit mm | Scale 1 : 1


Type	Machine			Ratchet		
Code	UNFDC 3006	UNFDC 3012	UNFDC 3018	UNFDW 3006	UNFDC 3012	UNFDC 3018

















Screw Driver

K Hex 1.2		Unit mm Scale 1 : 1	
Type	Machine		Ratchet
Code	MD 1230		HD 1220A
			

Solid Abutment Driver

· Only use UFil Narrow Kit


K D Ø4.0		● Sold separately	
		Unit mm Scale 1 : 1	
Length	6		12
Code	UNHDS 4006A		UNHDS 4012A
			

UFII Wide Kit

Kit Code I UF 09

- Kit composed of drill with outstanding cutting force and durability, etc.
- UF II Wide $\varnothing 5.9$ $\varnothing 6.4$ $\varnothing 6.9$ Fixture

Unit mm | Drill Scale 1 : 0.9



Surgical Tool

Trephine Drill



Scale 1 : 1.15

K				Unit mm Scale 1 : 1			
Code		TRD 4739		TRD 5244		TRD 5749	
							
Diameter 1		Ø3.9		Ø4.4		Ø4.9	
Diameter 2		Ø4.7		Ø5.2		Ø5.7	

Profile Drill



Scale 1 : 1.15

K				Unit mm Scale 1 : 1			
Fixture Size		Ø5.9		Ø6.4		Ø6.9	
Code		UPD 6010		UPD 6510		UPD 7010	
							
		Ø5.9		Ø6.4		Ø6.9	
Color		Yellow		Red		Green	

Surgical Tool

Final Drill



● Sold separately

Unit mm | Scale 1:1

K D Ø4.9						
Code	DSU 5007	DSU 5008	DSU 5010	DSU 5011	DSU 5013	DSU 5015
						
Length	7	8.5	10	11.5	13	15

K D Ø5.4						
Code	DSU 5507	DSU 5508	DSU 5510	DSU 5511	DSU 5513	DSU 5515
						
Length	7	8.5	10	11.5	13	15

K D Ø5.9						
Code	DSU 6007	DSU 6008	DSU 6010	DSU 6011	DSU 6013	DSU 6015
						
Length	7	8.5	10	11.5	13	15


Tap Drill

K


Unit mm | Scale 1 : 1

Fixture Size	Ø5.9	Ø6.4	Ø6.9
Code	UT(II) 6006	UT(II) 6506	UT(II) 7006

Color



Ø6.0



Ø6.5




Ø7.0

Fixture Driver


K D Ø3.6

Unit mm | Scale 1 : 1


Type	Machine			Ratchet		
Code	FDC 3506	FDC 3512	FDC 3518	FDW 3506A	FDW 3512A	FDW 3518A



6



12



18



6



12




18


100

UFII Surgical Protocol


Ø3.0 Fixture



Ø3.3 Fixture




Ø3.8 Fixture




UFII Surgical Protocol


Ø4.0 Fixture




Ø4.5 Fixture




Ø5.0 Fixture



Ø5.5 Fixture




Ø5.9 Fixture




UFII Surgical Protocol

Ø6.4 Fixture




Ø6.9 Fixture




UF Surgical Protocol

Ø3.8 Fixture




Ø4.0 Fixture




UF Surgical Protocol

Ø4.5 Fixture



Ø5.0 Fixture




DIOnavi. Kit

DIOnavi. Master Kit

Kit Code I UF 05

- Kit composed of drill with outstanding cutting force and durability, etc.
- UF II Narrow Ø3.0 Ø3.3 UF (II) Regular Ø3.8 Ø4.0 Ø4.5 Ø5.0 Fixture
- Guide Sleeve : GS 53

Unit mm I Drill Scale 1 : 0.9



Surgical Tool

Tissue Punch

- Gingiva can be removed stably using a fixed blade inside the punch.
(Flapless Surgery)



Unit mm | Scale 1 : 1


K D Ø3.0

Code

GTP 5125

Bone Flattening Drill

- Remaining gingiva is removed for easier drilling,
and irregular alveolar crest is flattened.



Unit mm | Scale 1 : 1

K D Ø5.1

Code

UBFD 5127

Guide Drill

- A hole on the bone is initially formed
on the bone for easier insertion of the initial drill.



Unit mm | Scale 1 : 1

K D Ø2.0 Length 5

Code

UGD 2005R

Drill Extension



Unit mm | Scale 1 : 1

K Length 11

Code

DE 3811

Surgical Tool

Drill Tube

· To fix the guide drill and the initial drill without shakage.



● Sold separately


Unit mm | Scale 1 : 1

Code	UDT 20		UDTE 20		UDTL 20	



Initial Drill

· A placement hole(osteotomy site) is formed on the cortical bone.



● Sold separately


Unit mm | Scale 1 : 1

Code	ISD 2005 ISD 2007 ISD 2008 ISD 2010 ISD 2011 ISD 2013 ISD 2015						




Profile Drill

- The cortical bone is further drilled in D1 and D2 bones, and excessive placement torque is prevented.

**K**

Unit mm | Scale 1 : 1

Fixture Size	Ø3.0	Ø3.3	Ø3.8	Ø4.0	Ø4.5	Ø5.0
Code	GPD 3005	GPD 3305	GPD 3805	GPD 4005	GPD 4505	GPD 5005



Bone Profile Drill


- When connecting abutment (healing) after implant placement, It is necessary to contour the bone for proper seating of the abutments. (to contour bone to allow for a wider emergence profile restoration)
- To flatten and remove bone that has grown around or over an implant which impedes its restoration. (When use a Abutment (with Healing Abutment))
- When abutment is completely tightened transmucosal recovery time is reduced.

● Sold separately

K

Unit mm | Scale 1 : 1

Code	BPD 4616GN	BPD 5608G
------	------------	-----------



Narrow




Regular

Type

Surgical Tool

Final Drill

· The drill hole is expanded to the final drilling.



K D Ø2.6							Sold separately
Unit mm Scale 1:1							
Code	USD 2705	USD 2707	USD 2708	USD 2710	USD 2711	USD 2713	USD 2715
							
Length	5	7	8.5	10	11.5	13	15

K D Ø2.9							Unit mm Scale 1:1
Code	USD 3005	USD 3007	USD 3008	USD 3010	USD 3011	USD 3013	USD 3015
							
Length	5	7	8.5	10	11.5	13	15

● Sold separately

K D Ø3.1

Unit mm | Scale 1 : 1

Code	USD 3205	USD 3207	USD 3208	USD 3210	USD 3211	USD 3213	USD 3215
							
Length	5	7	8.5	10	11.5	13	15

K D Ø3.7

Unit mm | Scale 1 : 1

Code	USD 3805	USD 3807	USD 3808	USD 3810	USD 3811	USD 3813	USD 3815
							
Length	5	7	8.5	10	11.5	13	15

K D Ø4.2

Unit mm | Scale 1 : 1

Code	USD 4305	USD 4307	USD 4308	USD 4310	USD 4311	USD 4313	USD 4315
							
Length	5	7	8.5	10	11.5	13	15

Surgical Tool


Tap Drill

- Excessive implant placement torque is prevented.
- $\varnothing 3.0/\varnothing 3.3$: Uses 5.3 Sleeve (GS 53BL: Blue)

● Sold separately

Unit mm | Scale 1 : 1

K	Fixture Size	$\varnothing 3.0$	$\varnothing 3.3$	$\varnothing 3.8$	$\varnothing 4.0$	$\varnothing 4.5$	$\varnothing 5.0$
	Code	GNTD 3015	GNTD 3315	GTD 3815	GTD 4015	GTD 4515	GTD 5015

						
	$\varnothing 3.0$	$\varnothing 3.3$	$\varnothing 3.8$	$\varnothing 4.0$	$\varnothing 4.5$	$\varnothing 5.0$
Color						

Abutment Profile Drill

- The cortical bone disrupting sitting of the abutment is removed, and emergence profile is formed.

K	Unit mm Scale 1 : 1		
Code	UAPD 5123	UAPD 5124	UAPD 5126

			
	$\varnothing 5.1$	$\varnothing 5.1$	$\varnothing 5.1$
Offset	9	10.5	12

Connector Extension

- Length of the connector is extended during implant placement.

K	Unit mm Scale 1 : 1			
Code	CE 6508A	CE 6508M	CE 6512A	CE 6512M

				
	12	12	12	12
Type	Ratchet	Machine	Ratchet	Machine




Implant Connector

- Implant is placed according to depth of the fixture and direction of the internal hex planned in advance.
- Recommended use of UF(II) Narrow : 20 times / Be careful about excessive torque of 50Ncm or above.

● Sold separately

Unit mm | Scale 1:1

K				Unit mm Scale 1 :
Type	Ratchet			Machine
Code	Non-Stopper	Stopper	Multi	Non-Stopper





UF II **Narrow**

GNIC 5309

GNIC 5309ST

GNIC 5317

GNICC 5309

UFII **Regular**

GIC 5309


GIC 5309ST

GIC 5317



GICC 5309

Surgical Tool

Screw Driver

K Hex 1.2	Unit mm Scale 1 : 1		
Code	MD 1222	MD 1230	HD 1215A
			
Type	Machine	Machine	Ratchet

Sqaure Screw Driver

K Hex 1.2	Unit mm Scale 1 : 1	
Code	SHD 1215	SHD 1220
		

Torque Wrench

· To place the fixture using an implant connector




K	
Code	DTW 0080

DIOnavi. Narrow Kit

Kit Code I UF 14

- Kit composed of drill with outstanding cutting force and durability, etc.
- UF(II) Narrow $\varnothing 3.0$ $\varnothing 3.3$ Fixture
- Guide Sleeve : GSL 36, GS36

Unit mm | Drill Scale 1 : 0.9



Surgical Tool

Tissue Punch

- Gingiva in the position of implant placement is incised into a round shape, and the inner double blade is used to remove the gum tissue cut into round shape.




Double blade
Gingiva can be removed stably using a fixed double blade inside the punch (flapless surgery).

Unit mm | Scale 1 : 1.15

K D Ø2.8	
Code	GTP 3430

Bone Flattening Drill

- Alveolar crest surface is flattened, and soft tissue remaining in the alveolar crest is removed after using the tissue punch.



Quadruple blade
A triple blade is used for stable drilling of the hard alveolar ridge.

Unit mm | Scale 1 : 1.15

K D Ø3.5	
Code	UBFD 3430

Guide Drill

- A hole in the bone is initially formed on the bone for easier insertion of the initial drill.




Double stepped drill
Stable drilling can be done due to short double step drill. This is a basic drill to maintain the accuracy of the path.

Unit mm | Scale 1 : 1.15

K D Ø2.0	
Code	UGD 2005N

Initial Drill

- It secures accurate depth and direction of drilling performed after achieving accurate hole formation during the initial steps.



K D Ø2.0						
Code	ISD 2007N	ISD 2008N	ISD 2010N	ISD 2011N	ISD 2013N	ISD 2015N

Drill Tube

- To fix the guide drill and the initial drill without shakage.



K D Ø3.5	
Code	UDTE 10

Profile Drill

- Used to further drill cortical bone in the mandible with bone density of D1 or D2 for the purpose of preventing excessive torque.

K			Unit mm Scale 1 : 1
Fixture Size	Ø3.0	Ø3.3	
Code	GNPD 3005DN	GNPD 3305DN	
Color			


Surgical Tool

Final Drill

● Sold separately


K D Ø2.5

Unit mm | Scale 1 : 1

Code	USD 2507DN	USD 2508DN	USD 2510DN	USD 2511DN	USD 2513DN	USD 2515DN
						
Length	7	8.5	10	11.5	13	15


K D Ø2.6

Unit mm | Scale 1 : 1

Code	USD 2707DN	USD 2708DN	USD 2710DN	USD 2711DN	USD 2713DN	USD 2715DN
						

K D Ø2.9


Unit mm | Scale 1 : 1

Code	USD 3007DN	USD 3008DN	USD 3010DN	USD 3011DN	USD 3013DN	USD 3015DN
						

Abutment Profile Drill

- Shape of abutment profile is formed by rotating the functional drill 360 degrees which removes the disruptive alveolar bone around the abutment or H-Scanbody.

K D Ø3.4	
Code	UAPD 3428




Unit mm | Scale 1 : 1.15

Tap Drill

- Uses 3.6 Sleeve (GSL 36YE : Yellow)

K Hex 1.7		Sold separately
Code	GNTD 3018DN	GNTD 3318DN

Unit mm | Scale 1 : 1



Color


Implant Connector

- If depth of placement and position of fixture hexagon are aligned with the guide sleeve projection during placement of the fixture, it is aligned with hex direction of the premade customized abutment.


K Hex 1.7		Unit mm Scale 1 : 1
Code	GNIC 3615	GNIC 3618



Type



Short



Long

Connector Extension

- Length of the connector is extended during implant placement.

K					 Sold separately
Unit mm Scale 1:1.2					
Code	CE 6508A	CE 6508M	CE 6512A	CE 6512M	

Unit mm | Scale 1 : 1.2



Ratchet



Machine



Ratchet



Machine




Type

Surgical Tool







Screw Driver Machine

· Used to connect the abutment screw (for handpiece) using an abutment screw driver for ratchet.






● Sold separately

K Hex 1.2		Unit mm Scale 1:1			
Code	MD 1219	MD 1222	MD 1230	MD 1234	MD 1239
     					

Screw Driver Ratchet

K Hex 1.2		Unit mm Scale 1:1				
Code	HD 1205A	HD 1210A	HD 1215A	HD 1220A	HD 1225A	HD 1230A
      						

Square Screw Driver


K Hex 1.2		Unit mm Scale 1:1			
Code	SHD 1210	SHD 1215	SHD 1220	SHD 1225	SHD 1230
    					

DIOnavi. Wide Kit

Kit Code | UF 11

- Kit composed of drill with outstanding cutting force and durability, etc.
- UFⅡ Regular Ø5.5 / UFⅡ Wide Ø5.9 Ø6.4 Fixture
- Guide Sleeve : GS 68GN

Unit mm | Drill Scale 1 : 0.9



Surgical Tool

Tissue Punch

· Gingiva can be removed stably using a fixed blade inside the punch

K

D Ø4.5

Code

GTP 6623



Unit mm | Scale 1 : 1.15

Bone Flattening Drill

· Irregular alveolar crest surface is flattened.

K

Code

UBFD 6623



Unit mm | Scale 1 : 1.15

Drill Tube

· Used for drilling of the guide drill and initial drill without oscillation

K

Code

UDT 53



Unit mm | Scale 1 : 1.2

Guide Drill

· A hole is initially formed on the bone for easier insertion of the initial drill.

K

D Ø2.0

Code

UGD 2002WUGD 2005W


Unit mm | Scale 1 : 1



Surgical Tool

Final Drill


· The drill hole is expanded to the final drilling.



K D Ø2.6					Unit mm Scale 1:1
Code	USD 2705DW	USD 2707DW	USD 2708DW	USD 2710DW	



K Length 5.5				Unit mm Scale 1:1
Code	USD 3205DW	USD 3805DW	USD 4305DW	



Surgical Tool

K D Ø4.7					Unit mm Scale 1 : 1	
Code	USD 4805DW	USD 4807DW	USD 4808DW	USD 4810DW		
						
Length	5.5	7	8.5	10		

K D Ø4.9					Unit mm Scale 1 : 1	
Code	USD 5005DW	USD 5007DW	USD 5008DW	USD 5010DW		
						
Length	5.5	7	8.5	10		

K D Ø5.4					Unit mm Scale 1 : 1	
Code	USD 5505DW	USD 5507DW	USD 5508DW	USD 5510DW		
						
Length	5.5	7	8.5	10		




Profile Drill

- Used to further drill cortical bone in the mandible with bone density of D1 or D2 for the purpose of preventing excessive torque

K		Unit mm Scale 1:1		
Fixture Size		Ø5.5	Ø5.9	Ø6.4
Code		GPD 5505DW	GPD 6005DW	GPD 6505DW
				
Color				

Tap Drill

- Excessive implant placement torque is prevented.

K		Unit mm Scale 1:1		
Fixture Size		Ø5.5	Ø5.9	Ø6.4
Code		GTD 5515DW	GTD 6015DW	GTD 6515DW
				
Color				

Connector Extension

- Length of the connector is extended during implant placement.

K		Unit mm Scale 1:1.1	
Code		CE 6508A	CE 6508M
			
Type		Ratchet	Machine

Abutment Profile Drill




- The cortical bone disrupting sitting of the abutment is removed, and emergence profile is formed.

K D Ø5.8		Unit mm Scale 1 : 1	
Code	UAPD 5822	UAPD 5823	
			
Offset	8	9.5	

Implant Connector






- Implant is placed according to depth of the fixture and direction of the internal hex planned in advance.
- Maximum Tightening Torque : 50Ncm

● Sold separately

K		Unit mm Scale 1 : 1		
Code	GWIC 6807ST	GWIC 6808ST	GWIC 6808	GWICC 6808
				
Type	Stopper	Stopper	Non-Stopper	Machine Non-Stopper
Offset	7	8	8 / 9.5	8 / 9.5

Square Screw Driver


● Sold separately

K Hex 1.2		Unit mm Scale 1 : 1			
Code	SHD 1210	SHD 1215	SHD 1220	SHD 1225	SHD 1230
					
	10	15	20	25	30

DIOnavi. Flapless Crestal Sinus Kit

Kit Code | SMK 02

Unit mm | Drill Scale 1 : 0.9



Surgical Tool




Ø2.0 Drill Stopper

- This is used to connect to the Ø2.0 initial drill of DIOnavi. Master to maintain constant depth during surgery.
- Range of length is 3~6mm with 1mm interval and anodizing and laser marking are done for each length.

Unit mm Scale 1 : 1					
K	Code	GST 2030BL	GST 2040YE	GST 2050RD	GST 2060BL
					

Sinus Round Drill

- Dome shaped tip is designed to prevent damaging of the membrane during drilling.
- Depth must be adjusted using a stopper.
- Non-irrigation and low-speed drilling (50rpm)

K D Ø3.2						
Code	GSRD 3217		GSRD 3218		GSRD 3219	GSRD 3221
						

Ø3.2 Drill Stopper For crestal approach

- This is used by connecting to the Ø3.2 round drill of Flapless Crestal Sinus Kit.
- Range of length is 3~13mm with 1mm interval to maintain constant depth during surgery.
- Anodizing and laser marking are done for each length.

Unit mm Scale 1 : 1							
K	Code	GST 03BL	GST 04YE	GST 05RD	GST 06BL	GST 07YE	GST 08RD
							
	Code	GST 09BL	GST 10YE	GST 11RD	GST 12BL	GST 13YE	
							

Depth Gauge & Bone Condenser

- After checking the thickness of the residual bone and whether the membrane is lifted, push in the bone into the lifted sinus.
- After separating the Guide, make sure to use it attaching the Stopper.

K	Unit mm
Code	GSDB 2619



Membrane Lifter

- Single use only

K	Unit mm
Code	WML 6627



Square Wrench

K	Unit mm
Code	DSW 050



Surgical Tool

Saline Tube

· Syringe ※ Separately sold product

K			● Sold separately Unit mm
	Code	MLT 40300	




DIONavi. Special Kit

Kit Code | SGF 02

· Kit composed of drill with outstanding cutting force and durability, etc.


Unit mm | Drill Scale 1 : 0.9



Surgical Tool

Anchor Drill

- A hole is formed to insert an anchor pin into the bone for the purpose of fixating the surgical guide.



K D Ø2.0	Unit mm Scale 1 : 1.15	
Code	AD 2019	

Stopper

- This is used when connecting to the Ø2.0 initial drill. Constant depth is maintained during surgery.





K D Ø2.9	Unit mm Scale 1 : 1.15	
Code	GST 2040YE	

Anchor Screw

- Used to connect the anchor screw using an anchor screw driver for ratchet

K D Ø2.0	Unit mm Scale 1 : 1	
Code	ASC 1515S	ASC 1519S



Anchor Screw Driver

- Used to connect the anchor screw using an anchor screw driver for ratchet



K	Unit mm Scale 1 : 1.15	
Code	ASD 2513	

Surgical Tool

Initial Drill



K D Ø2.0

Unit mm | Scale 1 : 1.15

Code

ISD 2015

Drill Tube



K

Unit mm | Scale 1 : 1.2

Code

UDT 25

Point Straight Drill

- Use of a drill tube (UDT 25) exclusive for the instrument (with a shape that prevents slipping on the bone) to form the guide hole in accurate position


K D Ø2.5

Unit mm | Scale 1 : 1

Code

PSD 2518


PSD 2525



Surgical Tool

Straight Drill

K Length 15		Unit mm Scale 1 : 1			
Code	USD 2715	USD 3015	USD 3215	USD 3815	USD 4315
Color					



Surgical Guide Fix Pin

- Used after Ø2.0 drilling to insert a surgical guide fix pin into the hole to fix the surgical guide in place and prevent shaking (side shaking is prevented after fixing)

K	Unit mm Scale 1 : 1.15	
Code	SGFP 2008	




Surgical Guide Fix

- Used to fix the surgical guide in place and prevent movement by connecting to the fixture placed earlier (side shaking is prevented after fixing)

K	Unit mm Scale 1 : 1		
Code	SGF 6309	SGF 6310	SGF 6312
Offset	 0	 1.5	 3
	9	10.5	12

Metal Needle

- Place the Metal needle tip deep into the osteotomy site, in order to prevent any damage caused by over-heating of the bone and complete removal of remaining tissues and particles at drilled bone cavity.
- Recommended capacity of the syringe : 30~50cc (Refrigerate before surgery)



O

Place the metal needle tip deep into the osteotomy site.



X

With the needle tip not all the way down to osteotomy it is hard to irrigate

K

Unit mm | Scale 1 : 1.15

Code	MNTE	MNTL
Height 1	30	25
Height 2	50	25

DIONavi. Drill Plate

- To increase efficiency and convenience of sterilizing drills and restoring them into the KIT set after DIONavi. Surgery
 - Sterilizable material
 - Size : 200 × 140 × 30 (mm)
- ※ Non-medical devices



K

Unit mm

Code	NP 01
------	-------

Surgical Kit

Sinus Master Kit

Kit Code | SMK 01


· Crestal Approach Technic

Sinus Technic based on low-speed drilling (50~100rpm) and stopper that perforates the inferior cortical wall and lifts the membrane safely and quickly

· Lateral Approach Technic

Round drill or core drill can be selected. The membrane can be lifted safely and quickly based on low-speed drilling (50~100rpm) and stopper.

Unit mm | Drill Scale 1 : 0.9




SURGICAL KIT

Sinus Master Kit


Surgical Tool Crestal

Ø2.0 Intial Drill



K D Ø2.0	Unit mm Scale 1 : 1.2	
Code	SID 2013	

Round Drill



K Length 13	Unit mm Scale 1 : 1.2	
Code	SRD 2813	SRD 3613

Stopper

· Stopper must be used for adjusting depth of the drill. It is distinguished by color and number.

K	Unit mm Scale 1 : 1										
Code	ST 02RD	03BL	04YE	05RD	06BL	07YE	08RD	09BL	10YE	11RD	12BL




Drilling Length	2	3	4	5	6	7	8	9	10	11	12
Stopper Length	15	14	13	12	11	10	9	8	7	6	5

Membrane Lifter

· Single use only



K	Unit mm Scale 1 : 1.2	
	Code	WML 6627



Depth Gauge

K	Unit mm	
	Code	SDG 2533




Bone Condenser

K	Unit mm	
	Code	SBC 2533




Surgical Tool Lateral

Round Drill

K		Unit mm Scale 1 : 1	
Code	SLRD 5504		SLRD 7004
 			

Core Drill

K		Unit mm Scale 1 : 1	
Code	SLCD 5504		SLCD 7004
 			

Stopper

· Stopper must be used for adjusting depth of the drill. It is distinguished by color and number.

K							Unit mm Scale 1 : 1	
Code	LST 05YE	10RD	15BL	20YE	25RD	30BL		
								
Drilling Length	0.5	1	2	1.5	2.5	3		
Stopper Length	7.5	7	6	6.5	5.5	5		

Surgical Protocol

Lateral Approach Technic


- Round drill or core drill can be selected. The membrane can be lifted and based on low-speed drilling (50~100rpm) and stopper.

Drilling Round Drill

Rounded tip to perform drilling

Caution

The stopper must be used. Irrigation at 1,200~1,400rpm




Drilling Core Drill

The edge of the rounded tip only can perform drilling
Beware of the remaining bone chips returning to original position after lifting of the maxillary sinus.


Caution

The stopper must be used. Irrigation at 600~800rpm



Lifting of the sinus membrane Sinus Lift

Lifting of the membrane on the lateral side



Surgical Protocol

Crestal Approach Technic

- Sinus Technic based on low-speed drilling (50~100rpm) and stopper that perforates the inferior cortical wall and lifts the membrane.


Drilling Initial Drill

Drilling is done 1~2mm short from the floor of the maxillary sinus measured by CT after attaching a stopper.

Caution

The stopper must be used to adjust depth.

Non-irritation and low-speed drilling at 50rpm




Sinus drill to perforate the bone underneath the maxillary sinus

Drilling is done to the same length as height of the residual bone of the maxillary sinus measured on CT after attaching a stopper.

Caution

The stopper must be used to adjust depth.

Non-irritation and low-speed drilling at 50rpm



Membrane lifter to lift the sinus membrane

The membrane lifter is used to inject saline solution into the osteotomy site.

Caution

Amount of injection measurement starts when pressure is felt, excluding the first 0.2cc~0.5cc before pressure is applied.

※ Amount of injection before pressure can differ according to height and expansion of the bone. For 1mm of sinus lifted, 0.1cc is injected.

Caution

If inferior border (A) of the maxillary sinus is well open

If pressure can be felt during injection of saline solution, the membrane is lifted and the pressure drops and saline solution injected.


If inferior border (A) of the maxillary sinus is not open well

After pressure felt from injection of saline solution, no more pressure can be applied or the nozzle is pushed out.

→ Drill the sinus drill 1mm deeper and retry.

Aspiration of saline solution with the nozzle kept in the hole


If negative pressure results after aspiration of injected saline solution mixed with blood, the membrane is fully lifted.



Sinus drill (2nd) to open and expand inferior border of the maxillary sinus.

After lifting the maxillary sinus membrane, the sinus drill is used to drilled 1mm deeper to expand the entrance of inferior border of the maxillary sinus.

Caution The stopper must be used to adjust depth.



Bone condenser to inject bone graft material

The bone condenser is used to push bone graft material into the maxillary sinus through the osteotomy site.

Caution

The stopper must be used to adjust depth.

Non-irrigation, low-speed drilling at 50rpm.

Caution


Determine volume of bone graft material.

When pressure is felt during injection of saline solution, the membrane is lifted.

Lifting height of the sinus membrane (mm)		1	2	3	4	5	6	7	8	9	10
Bone graft GBR (cc)	If implant is placed	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9	1.0
	If implant is not placed	0.3	0.6	0.9	1.2	1.5	1.8	2.1	2.4	2.7	3.0

Dispersion of bone graft material

Insert a depth gauge into the maxillary sinus and rotate it to evenly disperse bone graft material.



Final Drilling Final Drill

Drill 2mm deeper than depth of the sinus drill used earlier.

Caution Use a drill shorter than implant.



Placement of Implant

Implant that enters into the maxillary sinus pushes bone graft material away to disperse it. If amount of residual bone is 4mm or larger, satisfactory initial fixing force can be obtained. Also, temporary prosthesis can be restored immediately upon placement of implant.

Caution

If remaining bone is very thin with thickness of 3mm or less and initial fixing of implant fails, only bone graft on the maxillary sinus is performed without implant placement.




Dr.SOS+ Kit

Kit Code | DRS 02

- Problems like fracture of screw from excessive torque, wearing of connections, and jamming of foreign matters during prosthetic part of implant installation are resolved safely and quickly.

Unit mm | Drill Scale 1 : 0.9



Surgical Tool

Internal Guide

· Guide for accurate entry of the Remove Bar and Reverse Cutting Drill into the fixture

K UFII Narrow		Unit mm Scale 1:1	
Code		NINT 11H	NINT 11N
Type			
		Hex	Non-Hex

K UFII Regular / Wide			
Code		INT 11H	INT 11N
Type			
		Hex	Non-Hex

K Non-Sub. 8°			
Code		INT 8Octa	INT 8N
Type			
		Octa	Non-Octa


External Guide

K Hex		Unit mm Scale 1:1		
Code		EXT 24	EXT 27	EXT 34
				
				

Surgical Tool

Reverse Cutting Drill

- After inserting the end blade of the cutting drill in the fractured screw remains, remove the remains in a reverse rotation by friction.
- When it is difficult to remove with a reverse cutting drill, remove the screw remains by inserting a 1.2 hex removing driver on the hole formed on the screw remains and drilling in a reverse rotation




※Drill Speed : 1,200RPM / Reverse rotation
※In order to prevent high heat generation, osteotomy irrigation is necessary.

K		Unit mm Scale 1 : 1	
Code		RCD 32	RCD 34
Type		M1.4 / 1.6	M2.0 / 2.5

Remove Bar

- By attaching the end blade of the bar to the cross section of the fractured screw remains, remove the remains in a reverse rotation by friction



※ Drill Speed : 15rpm / Reverse rotation

K		Unit mm Scale 1 : 1	
Code		RB 32	RB 34
Type		M1.4 / 1.6	M2.0 / 2.5

Fixture Removing Driver

- Allows the removal of placed implant by inserting to the fixture and turning anti-clockwise.



Unit mm | Scale 1 : 1

Fixture Size	UFII Narrow / External NF		UFII Regular / Wide	
Code	NFRD 11S	NFRD 11L	FRD 11S	FRD 11L



Type

Short

Long

Short


Long

SURGICAL KIT

Dr.SOS+ Kit

1.2 Hex Removing Driver

- When use of the 1.2 Hex driver is not possible due to microorganisms in the 1.2 Hex connection part or the screw is stripped, use the removing driver for the removal of abutments or screws by inserting and turning anti-clockwise.




Unit mm | Scale 1 : 1

Code	HR 12S	HR 12L
------	--------	--------



0.8 Slot Driver

- If removal with the driver is not possible due to severe wear in the connection part, use the bur to create a slot with a minimal depth of 1mm and remove using the 0.8 Slot driver.




K	Unit mm Scale 1 : 1	
Code	SD 08	



Abutment Connection Removing Driver

- Used to remove a fractured abutment. Connect the screwdriver into the fixture and turn it counterclockwise.



K	Unit mm Scale 1 : 1	
Abutment Size	Narrow	Regular
Code	NHARD 2326	HARD 2326




Color



Screw Guide Tap

- Allows the restoral of distorted or damaged screw threads which hinder abutment-to-fixture attachment, caused by the penetration of microorganisms in the connection part or stripping of the screw thread.



K UF II Narrow / External Narrow

Unit mm | Scale 1 : 1

Code	SGT 1418A	SGT 1618A	SGT 2018A	SGT 2518A
				
Color				

※ This catalog contains all products of DIO Corporation.

Please check with your local sales office and your sales representative if you are able to order products.



Manufacturer

DIO Corporation

66, Centumseo-ro,
Haeundae-gu,
Busan, 48058, Republic of Korea